

Zircos-Com (Indirect)

Light cured, radiopaque nano-hybrid flowable composite material for the restorative treatment of anterior and posterior teeth.

Including Zirconia Filler



Application Ranges

- Class I-V fillings
- Reconstruction of traumatically affected anterior
- Correction of shape and shade for better aesthetic appearance
- locking, splinting of loose anterior
- Repairing veneer
- Filling deciduous teeth
- Core build-up under crowns
- Composite Inlays
- Depth of cure : 2.5mm (10 seconds)


Contents

- Kit-Zircos-Com 4g syr. *4ea
- A2, A3, B2, B3, Universal enamel E1, Universal enamel E2
- Zircos-Glaze Liqued


Refill Zircos-Com

- A1, A2, A3, A3.5, A4, A2E, A3E
- B1, B2, B3, C2, C3
- D2, D3, OA2, OA3, TL, WT, OP
- Universal enamel E1, Universal enamel E2


Physical properties



Compressive strength



Flexural strength



Elongation percentage